

SNF5 antibody

Cat. No. GTX33509

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:200
IP	Assay dependent

Not tested in other applications.

Calculated MW 44 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-376 of human SMARCB1 (NP_001007469.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX33509 IHC-P Image

IHC-P analysis of Human epithelioid sarcoma (ini-1 deletion) tissue using GTX33509 SNF5 antibody.

Antigen retrieval : 10 mM citrate buffer pH 6.0


Dilution : 1:200


GTX33509 IHC-P Image

IHC-P analysis of undifferentiated carcinoma of human esophagus tissue using GTX33509 SNF5 antibody.

Antigen retrieval : 10 mM citrate buffer pH 6.0

Dilution : 1:200


GTX33509 WB IP Image

IP analysis of HeLa cell lysates using GTX33509 SNF5 antibody.

Antibody amount : 3 µg/ 300 µg total extraction

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).