

p90 RSK1 (phospho Thr359/Ser363) antibody [RM233]

Cat. No. GTX33600

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000 - 1:2000
IHC-P	1:200 - 1:500

Not tested in other applications.

Calculated MW 83 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A phospho-peptide corresponding to human phospho-Rsk1 (Thr359/Ser363)
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX33600 WB Image

WB analysis of A431 cells treated (+) or nontreated (-) with EGF lysates using GTX33600 p90 RSK1 (phospho Thr359/Ser363) antibody [RM233].

Dilution : 1:1000



GTX33600 IHC-P Image

IHC-P analysis of human cerebral cortex tissue using GTX33600 p90 RSK1 (phospho Thr359/Ser363) antibody [RM233].

Dilution : 1:200



For full product information, images and publications, please visit our [website](#).