

EGFP antibody [G3]

Cat. No. GTX33910

Host	Human
Clonality	Monoclonal
Isotype	IgG2
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Species independent

References (2)

Package

100 µg

Applications

Application Note

Immunocytochemistry: 1:1000, Immunoprecipitation: 1:200

Product Note Recognizes EGFP, Enhanced Cyan Fluorescent Protein (ECFP) and Enhanced Yellow Fluorescent Protein (EYFP).

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Enhanced Green Fluorescent Proteins (EGFP).
Purification	Protein A-affinity purified.
Purity	≥95% by SDS-PAGE.
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).