

BAFF Receptor antibody [9B9]

Cat. No. GTX33920

Host	Rat
Clonality	Monoclonal
Isotype	IgG2b
Applications	FCM, Depletion
Reactivity	Mouse

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
Depletion	Assay dependent

Note : Depletion of B cells in vivo.

Not tested in other applications.

Product Note Recognizes mouse BAFF Receptor.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Y3 cells transfected with mouse BAFF-R.
Purification	Protein G-affinity purified.
Purity	≥95% by SDS-PAGE.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).