

Giantin antibody [TA10]

Cat. No. GTX33924

Host	Human
Clonality	Monoclonal
Isotype	IgG2
Applications	ICC/IF
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

Immunocytochemistry: 1:1000

Product Note Recognizes human and mouse giantin.

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human recombinant giantin.
Purification	Protein A-affinity purified.
Purity	??95% by SDS-PAGE.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX33924 Image



For full product information, images and publications, please visit our [website](#).