

LTBR antibody [3C8] (Azide free)

Cat. No. GTX33935-00

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	Activation
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Activation	Assay dependent
------------	-----------------

Note : Agonist inducing BAFF, chemokines and integrins in vitro and in vivo.

Not tested in other applications.

Product Note Recognizes mouse LTBR.

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Immunogen Recombinant mouse LTBR (cysteine-rich region/aa 31-221).

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).