

## RAB6 (Q72L mutant specific) antibody [AA2] (ATTO-488)

**Cat. No. GTX33942-13**

<b>Host</b>	Human
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2
<b>Applications</b>	ICC/IF
<b>Reactivity</b>	Human, Mouse, Drosophila

**Package**  
100 µg

## Applications

**Application Note**

Immunocytochemistry: 1:1000


**Product Note** Recognizes human, mouse and Drosophila GTP-bound Rab6a and Rab6b and mutant Rab6Q72L. Does not detect Rab6-GDP.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant Rab6AQ72L, a GTP-locked mutant of Rab6A in which Gln72 is replaced by Leu.
<b>Purification</b>	Protein A-affinity purified.
<b>Conjugation</b>	ATTO-488 <a href="#">Wavelength</a>

**Note**  
For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  
  
Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX33942-13 Image



For full product information, images and publications, please visit our [website](#).