

RAB6 (Q72L mutant specific) antibody [AA2] (ATTO-488)

Cat. No. GTX33942-13

Host	Human
Clonality	Monoclonal
Isotype	IgG2
Applications	ICC/IF
Reactivity	Human, Mouse, Drosophila

Package
100 µg

Applications

Application Note

Immunocytochemistry: 1:1000

Product Note

Recognizes human, mouse and Drosophila GTP-bound Rab6a and Rab6b and mutant Rab6Q72L. Does not detect Rab6-GDP.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Rab6AQ72L, a GTP-locked mutant of Rab6A in which Gln72 is replaced by Leu.
Purification	Protein A-affinity purified.
Conjugation	ATTO-488 Wavelength

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX33942-13 Image



For full product information, images and publications, please visit our [website](#).