

DR5 antibody [HS201] (Azide free)

Cat. No. GTX33945

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, FCM, IP, Neutralizing/Inhibition
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
FCM	Assay dependent
IP	Assay dependent
Neutralizing/Inhibition	Assay dependent

Note : Blocks TRAIL-R2 mediated killing if applied in solution)

Not tested in other applications.

Product Note Recognizes human TRAIL-R2. Does not cross-react with human TRAIL-R1, -R3 or -R4.

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human TRAIL-R2 (DR5).
Purification	Protein G-affinity purified.
Purity	??95% by SDS-PAGE.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).