

HSL antibody

Cat. No. GTX33956

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

WB: 1:500 - 1:1000


Calculated MW	117 kDa. (Note)
Product Note	This antibody detects endogenous levels of HSL protein.

Properties


Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 520-570 of Human HSL.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX33956 IHC-P Image

Immunohistochemistry (IHC) analyzes of HSL antibody in paraffin-embedded human colorectal cancer carcinoma tissue at 1:100.


GTX33956 WB Image

Western blot (WB) analysis of HSL antibody at 1:500 dilution Lane1:HEK293T whole cell lysate; Lane2:Raw264.7 whole cell lysate; Lane3:PC12 whole cell lysate



For full product information, images and publications, please visit our [website](https://www.genetex.com).