

PPAR delta antibody [2F9]

Cat. No. GTX34133

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 0.5% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant peptide derived from PPAR Delta
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX34133 IHC-P Image


IHC-P analysis of human brain tissue using GTX34133 PPAR delta antibody [2F9].
Dilution : 1:200


GTX34133 IHC-P Image


IHC-P analysis of human uterus tissue using GTX34133 PPAR delta antibody [2F9]. Negative control (the lower left coner) was secondary antibody only.
Antigen retrieval : Sodium citrate pH6.0 was used for antibody retrieval (>98°C, 20min)
Dilution : 1:200


GTX34133 IHC-P Image

IHC-P analysis of human lung tissue using GTX34133 PPAR delta antibody [2F9]. Negative control (the lower left coner) was secondary antibody only.
Antigen retrieval : Sodium citrate pH6.0 was used for antibody retrieval (>98°C, 20min)
Dilution : 1:200


GTX34133 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX34133 PPAR delta antibody [2F9].
Dilution : 1:200


GTX34133 IHC-P Image

IHC-P analysis of human liver tissue using GTX34133 PPAR delta antibody [2F9]. Negative control (the lower left coner) was secondary antibody only.
Antigen retrieval : Sodium citrate pH6.0 was used for antibody retrieval (>98°C, 20min)
Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).