

beta Catenin antibody [4F2]

Cat. No. GTX34339

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:2000
IHC-P	1:200-1:500

Not tested in other applications.

Calculated MW 86 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.5% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Protein
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX34339 WB Image


WB analysis of various sample lysates using GTX34339 beta Catenin antibody [4F2].

Lane 1 : HeLa cell lysate

Lane 2 : 293T cell lysate

Lane 3 : Mouse Liver tissue lysate

Lane 4 : Rat Liver tissue lysate




GTX34339 IHC-P Image

IHC-P analysis of human liver cancer tissue using GTX34339 beta Catenin antibody [4F2]. Negative control (the lower left coner) was secondary antibody only.


Antigen retrieval : Sodium citrate pH6.0 was used for antibody retrieval (>98°C, 20min)

Dilution : 1:200



GTX34339 IHC-P Image

IHC-P analysis of human lung carcinoma tissue using GTX34339 beta Catenin antibody [4F2].



GTX34339 IHC-P Image

IHC-P analysis of human liver tissue using GTX34339 beta Catenin antibody [4F2]. Negative control (the lower left coner) was secondary antibody only.

Antigen retrieval : Sodium citrate pH6.0 was used for antibody retrieval (>98°C, 20min)

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).