

M-CSF Receptor antibody

Cat. No. GTX34385

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

WB: 1:500-1:1000, IHC-P: 1:50-1:200

Calculated MW 108 kDa. ([Note](#))


Product Note c-Fms (Q718) pAb detects endogenous levels of c-Fms protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 691-740 of Human c-Fms.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX34385 IHC-P Image

Immunohistochemistry (IHC) analyzes of c-Fms pAb in paraffin-embedded human colorectal carcinoma tissue at 1:50.



GTX34385 WB Image

Western blot (WB) analysis of c-Fms pAb at 1:500 dilution Lane1:Hela cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate



For full product information, images and publications, please visit our [website](https://www.genetex.com).