

AMACR antibody [AMACR/1723]

Cat. No. GTX34410

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
IHC-P	1-2µg/ml for 30 minutes at RT

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 42 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human AMACR protein
Purification	Protein A/G purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX34410 WB Image


WB analysis of kidney tissue lysate using GTX34410 AMACR antibody [AMACR/1723].


GTX34410 IHC-P Image

IHC-P analysis of human prostate carcinoma tissue using GTX34410 AMACR antibody [AMACR/1723].


GTX34410 WB Image

WB analysis of PC3 cell lysate using GTX34410 AMACR antibody [AMACR/1723].


GTX34410 WB Image

WB analysis of liver tissue lysate using GTX34410 AMACR antibody [AMACR/1723].



For full product information, images and publications, please visit our [website](https://www.genetex.com).