

Chromogranin A antibody [LK2H10]

Cat. No. GTX34604

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	WB, IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat, Pig, Monkey

Reference (2)

Package

100 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	0.1-0.2µg/ml for 30 min at RT
IHC-Fr	Assay dependent

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW	51 kDa. (Note)
Product Note	We do not recommend use of this product for Rabbit, Sheep, Guinea pig samples.

PROPERTIES

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human pheochromocytoma
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2024 / 05 / 20 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX34604 IHC-P Image

IHC-P analysis of mouse pancreas tissue using GTX34604 Chromogranin A antibody [LK2H10].



GTX34604 IHC-P Image

IHC-P analysis of human small cell lung carcinoma tissue using GTX34604 Chromogranin A antibody [LK2H10].



GTX34604 IHC-P Image

IHC-P analysis of human adrenal gland tissue using GTX34604 Chromogranin A antibody [LK2H10].



For full product information, images and publications, please visit our [website](#).