

## Chromogranin A antibody [LK2H10]

**Cat. No. GTX34604**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Reactivity</b>	Human, Mouse, Rat, Pig, Monkey

References ( 2 )

Package

100 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	0.1-0.2µg/ml for 30 min at RT
IHC-Fr	Assay dependent

**Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.**

Not tested in other applications.

<b>Calculated MW</b>	51 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	We do not recommend use of this product for Rabbit, Sheep, Guinea pig samples.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.05% BSA
<b>Preservative</b>	0.05% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Human pheochromocytoma
<b>Purification</b>	Protein A/G purified
<b>Conjugation</b>	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


## DATA IMAGES

**GTX34604 IHC-P Image**

IHC-P analysis of mouse pancreas tissue using GTX34604 Chromogranin A antibody [LK2H10].

**GTX34604 IHC-P Image**

IHC-P analysis of human small cell lung carcinoma tissue using GTX34604 Chromogranin A antibody [LK2H10].

**GTX34604 IHC-P Image**

IHC-P analysis of human adrenal gland tissue using GTX34604 Chromogranin A antibody [LK2H10].



For full product information, images and publications, please visit our [website](#).