

Cytochrome C antibody [CTC05]

Cat. No. GTX34620

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2b
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat, Drosophila, Dog, Frog, Horse, pigeon

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2μg/ml
ICC/IF	1-2μg/ml
IHC-P	0.25-0.5μg/ml for 30 minutes at RT
FCM	1-2μg/10 ⁶ cells

Note: Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 12 kDa. (Note)

Properties	
Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant cytochrome c protein
Purification	Protein A/G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 2




For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX34620 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX34620 Cytochrome C antibody [CTC05].


Green: Primary antibody

Blue: nucleus



GTX34620 IHC-P Image

IHC-P analysis of human heart tissue using GTX34620 Cytochrome C antibody [CTC05].



GTX34620 FCM Image

FACS analysis of PFA fixed HeLa cells using GTX34620 Cytochrome C antibody [CTC05].

Blue: Primary antibody Red: Isotype control



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 2 of 2