

Cytokeratin 8 antibody [K8.8]

Cat. No. GTX34670

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
ICC/IF	1-2µg/ml
IHC-P	1-2µg/ml for 30 min at RT
FCM	1-2µg/10 ⁶ cells

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Product Note Epitope of this MAb is located between aa343-357.
We do not recommend use of this product for Rat samples.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Keratin preparation from a human carcinoma
Purification	Protein A/G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX34670 ICC/IF Image

ICC/IF analysis of HCT116 cells using GTX34670 Cytokeratin 8 antibody [K8.8].

Green : Primary antibody


Red : nucleus


GTX34670 FCM Image

FACS analysis of HeLa cells using GTX34670 Cytokeratin 8 antibody [K8.8].

Blue : Primary antibody

Red : Isotype control


GTX34670 WB Image

WB analysis of HCT116 cell lysate using GTX34670 Cytokeratin 8 antibody [K8.8].



For full product information, images and publications, please visit our [website](https://www.genetex.com).