

Fascin 1 antibody [FSCN1/417]

Cat. No. GTX34718

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
ICC/IF	1-2µg/ml
IHC-P	1-2µg/ml for 30 minutes at RT
FCM	1-2µg/10 ⁶ cells

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length recombinant human FSCN1 protein
Purification	Protein A/G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx34718 WB Image

WB analysis of HeLa cell lysate using GTx34718 Fascin 1 antibody [FSCN1/417].




GTx34718 FCM Image

FACS analysis of PFA-fixed K562 cells using GTx34718 Fascin 1 antibody [FSCN1/417].

Blue : Primary antibody

Red : Isotype control



GTx34718 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx34718 Fascin 1 antibody [FSCN1/417].

Green : Primary antibody

Red : Phalloidin

Blue : nuclear counterstain



For full product information, images and publications, please visit our [website](https://www.genetex.com).