

FOXA1 antibody [FOXA1/1515]

Cat. No. GTX34738

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
ICC/IF	0.5-1µg/ml
IHC-P	1-2µg/ml for 30 minutes at RT
FCM	0.5-1µg/10 ⁶ cells

Note : Staining of formalin-fixed tissues is enhanced by heating tissue sections in 10mM Tris buffer with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Product Note We do not recommend use of this product for Rat samples.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human FOXA1 protein fragment (aa372-472) (exact sequence is proprietary)
Purification	Protein A/G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX34738 WB Image

WB analysis of HeLa cell lysate using GTX34738 FOXA1 antibody [FOXA1/1515].


GTX34738 IHC-P Image


IHC-P analysis of human prostate carcinoma tissue using GTX34738 FOXA1 antibody [FOXA1/1515].


GTX34738 FCM Image

FACS analysis of PFA-fixed MCF-7 cells using GTX34738 FOXA1 antibody [FOXA1/1515].

Blue : Primary antibody

Red : Isotype control


GTX34738 IHC-P Image

IHC-P analysis of human prostate carcinoma tissue using GTX34738 FOXA1 antibody [FOXA1/1515].



For full product information, images and publications, please visit our [website](https://www.genetex.com).