

## MALT1 antibody [SPM578]

Cat. No. GTX34813

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG1       |
| Applications | WB, IHC-P  |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution          |
|--------------------|-------------------------------|
| WB                 | 1-2µg/ml                      |
| IHC-P              | 1-2µg/ml for 30 minutes at RT |

**Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.**

Not tested in other applications.

**Calculated MW** 92 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 0.05% BSA                                                                                                                                                                                                             |
| Preservative  | 0.05% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.2 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| Immunogen     | Human MALT1 recombinant fragment (aa701-808)                                                                                                                                                                               |
| Purification  | Protein A/G purified                                                                                                                                                                                                       |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX34813 IHC-P Image**

IHC-P analysis of human tonsil tissue using GTX34813 MALT1 antibody [SPM578].



For full product information, images and publications, please visit our [website](#).