

# MUC1 antibody [rMUC1/960]

**Cat. No. GTX34853**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Mouse      |
| <b>Clonality</b>    | Monoclonal |
| <b>Isotype</b>      | IgG2b      |
| <b>Applications</b> | IHC-P      |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                                                                                                                                                                 | Recommended dilution          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|
| IHC-P                                                                                                                                                                              | 1-2µg/ml for 30 minutes at RT |
| <b>Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.</b> |                               |

Not tested in other applications.

### Product Note

This MAb reacts with MUC1, a large transmembrane glycoprotein expressed on the ductal surface of normal glandular epithelia. The dominant epitope of this MAb involves both amino acids as well as sugar moieties. Neuraminidase treatment destroys the antigen. It is a very good tracer agent in CA15.3 assays. The extracellular domain of MUC1 largely consists of a highly conserved, O-glycosylated 20 amino acids tandem repeat which can occur 30-100 times per molecule depending on the length of the allele involved. In the vast majority of human carcinomas this protein is up-regulated and poorly glycosylated and appears on the cell surface in a non-polarized fashion.

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 0.05% BSA                                                                                                                                                                                                             |
| <b>Preservative</b>  | 0.05% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.2 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| <b>Immunogen</b>     | Recombinant human MUC1 protein                                                                                                                                                                                             |
| <b>Purification</b>  | Protein A/G purified                                                                                                                                                                                                       |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

#### DATA IMAGES



#### GTx34853 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTx34853 MUC1 antibody [rMUC1/960].



For full product information, images and publications, please visit our [website](#).