

Golden Syrian Hamster IgG isotype control (FITC)

Cat. No. GTX35006**Host** Golden Syrian Hamster**Clonality** Polyclonal**Isotype** IgG**Applications** FCM, ELISA**Package**

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution **Recommended dilution**

FCM Assay dependent

ELISA Assay dependent

Not tested in other applications.

Properties

Form Liquid**Buffer** PBS**Preservative** 0.1% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Conjugation** Fluorescein isothiocyanate (FITC) [Wavelength](#)**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).