

Golden Syrian Hamster IgG isotype control (PE)

Cat. No. GTX35008

Host	Golden Syrian Hamster
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA

Package 50 μg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Properties	
Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Phycoerythrin (PE) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 16 Page 1 of 1