

Mouse IgG isotype control (FITC)

Cat. No. GTX35010

Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA

References (1)

Package

50 µg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).