

## Mouse IgG2c isotype control [6.3] (Biotin)

Cat. No. GTX35013

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2c
Application	FACS, ELISA

Reference ( 2 )

Package

50 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

## PROPERTIES

Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Conjugation	Biotin
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).