

## Mouse IgG1 isotype control [15H6] (FITC)

Cat. No. GTX35015

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM, ELISA

References ( 8 )  
Package  
50 test

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10µl/10*6 cells
ELISA	Assay dependent

Not tested in other applications.

**Product Note** T-2 mycotoxin

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	T-2 mycotoxin
Conjugation	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).