

Mouse IgG2a isotype control [HOPC-1] (PE)

Cat. No. GTX35022

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM, ELISA

References (1)

Package

50 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10 μ l/10*6 cells
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Conjugation	Phycoerythrin (PE) Wavelength
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).