

Mouse IgG2b isotype control [A-1] (FITC)

Cat. No. GTX35024

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2b
Applications	FCM, ELISA, IHC

References (1) Package 50 test

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Chicken IgA

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Immunogen	Chicken IgA
Conjugation	Fluorescein isothiocyanate (FITC) <u>Wavelength</u>
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 04 Page 1 of 1