

# Mouse IgG2b isotype control [A-1] (FITC)

**Cat. No. GTX35024**

|                     |                 |
|---------------------|-----------------|
| <b>Host</b>         | Mouse           |
| <b>Clonality</b>    | Monoclonal      |
| <b>Isotype</b>      | IgG2b           |
| <b>Applications</b> | FCM, ELISA, IHC |

References ( 1 )  
 Package  
 50 test

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| FCM                | Assay dependent      |
| ELISA              | Assay dependent      |
| IHC                | Assay dependent      |

Not tested in other applications.

**Product Note**      Chicken IgA

## Properties

|                     |   |
|---------------------|---|
| <b>Form</b>         | Liquid  |
| <b>Buffer</b>       | PBS   |
| <b>Preservative</b> | 0.1% Sodium azide   |
| <b>Storage</b>      | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light. |
| <b>Immunogen</b>    | Chicken IgA   |
| <b>Conjugation</b>  | Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>  |

**Note**  
 For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  
 Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).