

Mouse IgG2b isotype control [A-1] (FITC)

Cat. No. GTX35024

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	FACS, ELISA, IHC

Reference (1)

Package

50 test

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Chicken IgA

PROPERTIES

Form Liquid

Buffer PBS

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.

Immunogen Chicken IgA

Conjugation Fluorescein isothiocyanate (FITC)

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).