

Mouse IgG3 isotype control [B10]

Cat. No. GTX35027

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Application	ICC/IF, FACS, IP, ELISA, in vitro

Reference (2)
Package
100 µg

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Unknown

PROPERTIES

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).