

Mouse IgG3 isotype control [B10] (FITC)

Cat. No. GTX35028

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Application	FACS, ELISA

Reference (2)
Package
50 test

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Unknown

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Conjugation	Fluorescein isothiocyanate (FITC)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).