

Mouse IgA isotype control [S107]

Cat. No. GTX35045

Host	Mouse
Clonality	Monoclonal
Isotype	IgA
Applications	FCM, ELISA

References (3)

Package

100 µg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Unknown

Properties

Form Liquid**Buffer** Borate buffered saline**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).