

Rat IgG isotype control

Cat. No. GTX35047

Host	Rat
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA

References (2)

Package

100 µg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	5 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).