

Rat IgG1 isotype control [KLH/G1-2-2] (PE)

Cat. No. GTX35054

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM, ELISA

References (1)

Package

50 µg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Keyhole limpet hemocyanin (KLH)

Properties

Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Keyhole limpet hemocyanin (KLH)
Conjugation	Phycoerythrin (PE) Wavelength
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).