

Rat IgG1 isotype control [KLH/G1-2-2] (PE)

Cat. No. GTX35054

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Application	FACS, ELISA

Reference (1)
Package
50 µg

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Keyhole limpet hemocyanin (KLH)

PROPERTIES

Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Keyhole limpet hemocyanin (KLH)
Conjugation	Phycoerythrin (PE)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).