

Rat IgG2a isotype control [KLH/G2a-1-1]

Cat. No. GTX35055

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Application	FACS, ELISA

Reference (1)
Package
100 µg

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Keyhole limpet hemocyanin (KLH)

PROPERTIES

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Keyhole limpet hemocyanin (KLH)
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).