

Rat IgG2c isotype control [SB68b] (Biotin)

Cat. No. GTX35064

Host	Rat
Clonality	Monoclonal
Isotype	IgG2c
Applications	FCM, ELISA

References (1)

Package

100 µg

Applications

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Keyhole limpet hemocyanin (KLH)

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Keyhole limpet hemocyanin (KLH)
Conjugation	Biotin
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).