

# Rat IgM isotype control [NIP/M-2] (PE)

**Cat. No. GTX35077**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgM
<b>Applications</b>	FCM, ELISA

**Package**

50 µg

## Applications

### Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

**Product Note** NIP

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, stabilizer
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	NIP
<b>Conjugation</b>	Phycoerythrin (PE) <a href="#">Wavelength</a>
<b>Note</b>	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).