

Mouse IgA isotype control [S107] (PE)

Cat. No. GTX35080

| | |
|--------------------|-------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgA |
| Application | FACS, ELISA |

Reference (3)

Package

50 test

APPLICATION

Application Note

Use at a concentration/dilution equivalent to that of your experimental primary antibody.

Product Note Unknown

PROPERTIES

Form Liquid

Buffer PBS, stabilizer

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.

Conjugation Phycoerythrin (PE)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).