

SUMO1 antibody [SM1/495]

Cat. No. GTX35092

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2µg/ml
ICC/IF	1-2µg/ml
IHC-P	1-2µg/ml for 30 minutes at RT
FCM	1-2µg/10 ⁶ cells

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW 12 kDa. ([Note](#))

Product Note This MAb is specific to SUMO-1 and shows no cross-reaction with either SUMO-2 or SUMO-3.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human SUMO1 protein
Purification	Protein A/G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX35092 ICC/IF Image

ICC/IF analysis of paraformaldehyde-fixed HepG2 cells using GTX35092 SUMO1 antibody [SM1/495].

Green : Primary antibody


Red : Phalloidin


GTX35092 FCM Image

FACS analysis of paraformaldehyde-fixed HepG2 cells using GTX35092 SUMO1 antibody [SM1/495].

Blue : Primary antibody

Red : Isotype control


GTX35092 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX35092 SUMO1 antibody [SM1/495].



For full product information, images and publications, please visit our [website](https://www.genetex.com).