

RHOBTB3 antibody

Cat. No. GTX35184

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Rat

Package
100 µl

Applications

Application Note

WB: 1:500-1:1000, IHC: 1:50-1:200


Calculated MW 69 kDa. ([Note](#))

Product Note RHOBTB3 (L365) polyclonal antibody detects endogenous levels of RHOBTB3 protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 341-390 of Human RHOBTB3.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX35184 WB Image

Western blot (WB) analysis of RHOBTB3 (L365) polyclonal antibody at 1:1000 dilution Lane1:HKE293T whole cell lysate(56ug) Lane2:MCF-7 whole cell lysate(48ug) Lane3:The Pancreas tissue lysate of Rat(43ug)



For full product information, images and publications, please visit our [website](https://www.genetex.com).