

BRG1 antibody

Cat. No. GTX35204

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:20 - 1:50
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 185 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide of human SMARCA4
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX35204 WB Image

WB analysis of normal (control) and knockout (KO) 293T cell lysate using GTX35204 BRG1 antibody.

Dilution : 1:500


Loading : 25µg per lane



GTX35204 IHC-P Image

IHC-P analysis of human colon tissue using GTX35204 BRG1 antibody.


Dilution : 1:100



GTX35204 IHC-P Image

IHC-P analysis of human colon carcinoma tissue using GTX35204 BRG1 antibody.

Dilution : 1:100



GTX35204 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX35204 BRG1 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).