

UBE2O antibody

Cat. No. GTX35205

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000 - 1:4000
ICC/IF	1:50 - 1:200
IP	1:50 - 1:200

Not tested in other applications.

Calculated MW 141 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1073-1292 of human UBE2O (NP_071349.3).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX35205 IP Image

IP analysis of A431 cell lysate using GTX35205 UBE2O antibody.

Antibody amount : 3 μ g / 200 μ g lysate


Dilution : 1:1000


GTX35205 WB Image

WB analysis of various sample lysates using GTX35205 UBE2O antibody.

Dilution : 1:1000

Loading : 25 μ g per lane


GTX35205 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX35205 UBE2O antibody.

Blue : DAPI

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).