

# MEK1 antibody

**Cat. No. GTX35212**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, IHC-P         |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000 - 1:3000      |
| IHC-P              | 1:50 - 1:200         |

Not tested in other applications.

**Calculated MW** 43 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAP2K1 (NP_002746.1).                                                                                                                    |
| <b>Purification</b>  | Purified by affinity chromatography                                                                                                                                                                                        |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |


### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX35212 WB Image**

WB analysis of various sample lysates using GTX35212 MEK1 antibody. The signal was developed with ECL plus-Enhanced.


Dilution : 1:3000

Loading : 25µg per lane


**GTX35212 IHC-P Image**


IHC-P analysis of rat brain tissue using GTX35212 MEK1 antibody.

Dilution : 1:100


**GTX35212 IHC-P Image**

IHC-P analysis of mouse kidney tissue using GTX35212 MEK1 antibody.

Dilution : 1:100


**GTX35212 WB Image**

WB analysis of normal (control) and knockout (KO) HeLa cell lysate using GTX35212 MEK1 antibody.

Dilution : 1:1000

Loading : 25µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).