

## IFITM3 antibody

Cat. No. GTX35214

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | 1:50 - 1:100         |
| IHC-P              | 1:50 - 1:200         |

Not tested in other applications.

Calculated MW 15 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Recombinant fusion protein containing a sequence corresponding to amino acids 1-133 of human IFITM3 (NP_066362.2).                                                                                                         |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX35214 ICC/IF Image

ICC/IF analysis of L929 cells using GTX35214 IFITM3 antibody.

Blue : DAPI

Dilution : 1:100




### GTX35214 ICC/IF Image

ICC/IF analysis of U2OS cells using GTX35214 IFITM3 antibody.

Blue : DAPI

Dilution : 1:100



### GTX35214 IHC-P Image

IHC-P analysis of human colon carcinoma tissue using GTX35214 IFITM3 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](#).