

BAF57 antibody

Cat. No. GTX35217

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IP	1:20 - 1:50

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-411 of human SMARCE1 (NP_003070.3).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX35217 WB Image


WB analysis of various sample lysates using GTX35217 BAF57 antibody.

Dilution : 1:1000

Loading : 25µg per lane


GTX35217 ICC/IF Image


ICC/IF analysis of HeLa cells using GTX35217 BAF57 antibody.


GTX35217 IP Image

IP analysis of Jurkat cell lysate using GTX35217 BAF57 antibody.

Antibody amount : 3µg / 150µg lysate

Dilution : 1:1000


GTX35217 WB Image

WB analysis of normal (control) and knockout (KO) HeLa cell lysate using GTX35217 BAF57 antibody.

Dilution : 1:1000

Loading : 25µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).