

CD163 antibody

Cat. No. GTX35247

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human, Mouse, Rat, Pig

References (10)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1000
IHC-P	1:50-400
FCM	Assay dependent

Not tested in other applications.

Calculated MW 125 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human CD163(800-850).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX35247 IHC-P Image

IHC-P analysis of human testis tissue using GTX35247 CD163 antibody.
Dilution : 1:100




GTX35247 IHC-P Image

IHC-P analysis of human testis tissue using GTX35247 CD163 antibody.
Dilution : 1:100



GTX35247 IHC-P Image

IHC-P analysis of mouse intestine tissue using GTX35247 CD163 antibody.
Dilution : 1:200



GTX35247 IHC-P Image

IHC-P analysis of rat colon tissue using GTX35247 CD163 antibody.



For full product information, images and publications, please visit our [website](#).