

CDC14B antibody

Cat. No. GTX36294

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1-4 µg/ml
IHC-P	1:20- 1:50

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 40% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	VENQDQQEPEPYSDDDEINGVTQGDRSRAL
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX36294 IHC-P Image

IHC-P analysis of human cerebellum tissue using GTX36294 CDC14B antibody. It shows strong nuclear membranous positivity in Purkinje cells.


GTX36294 ICC/IF Image

ICC/IF analysis of A431 cells using GTX36294 CDC14B antibody. It shows positivity in nuclear membrane.



For full product information, images and publications, please visit our [website](#).